



Total Intense & PowerGo

Technical Specifications

Product name:	Total Intense v2 & PowerGo v2
Material	Aluminium
Weight	Total Intense: 33 g (Incl. Framepart & Cable) PowerGo: 155 g
Brightness settings	Low Medium High
Operating time	16h 12.5h 9h
Charging time	Approx. 3.5-5 hours
Light intensity	At 25 cm: 51,500 lux At 35 cm: 25,400 lux
Light spot at working distance 40 cm	63 mm
Light colour	6,500 K
Colour rendering index	>70CRI
Cable type	Lamp Cable: Custom USB-C, 1.3 m Power Cord: USB-C 5W 2A
Battery type	21,700 Lithium-ion rechargeable battery
Battery rating	In: 5 VDC 2.0-2.5 A Out: 3 VDC 0.5-2.0 A
Warranty	Total Intense: 3 years PowerGo: 1 year
Product information	Flicker Free Light Rotatable belt clip Easy-to-operate Power button Over-current protection Over-voltage protection