



Kepler Kompakt

Technical Specifications

Variant	Kepler Kompakt v2 – 3.4x, 5.0x & 7.6x		
Lenses	Plano or Prescription		
BlueLight Protection (BLP)	Yes (optional add-on)		
Prisma Lens System	Yes (optional add-on)		
Fine-focus available	Yes		
Magnification	3.4x	5.0x	7.6x
Working distance	3.4x & 5.0x:	From 33 cm to 65 cm	
	7.6x:	From 39 cm to 65 cm	
Field of view (at 45 cm)	3.4x: 9.5 cm	5.0x: 6.2 cm	7.6x: 4.3 cm
Depth of field (up to)	3.4x: 11.5 cm	5.0x: 8.0 cm	7.6x: 3.5 cm
Mount System	TTL (Through the lens)		
Weight of two oculars with correction lenses	3.4x: 37.3 g	5.0x: 39.2 g	7.6x: 43.7 g
Loupe colour	Anthracite grey/Metallic grey		
Loupe material	Carbon fibre reinforced composite and Titanium		
Frames available	All		