



## Kepler Advanced Technical Specifications

Loupe lens system	Kepler
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Power-Up System	No, not needed (comes with 4 mags in one)
Magnifications	Turn'N'Click option of: 4.0x + 5.0x + 6.0x + 7.0x
Working distance	From 31 to 69 cm
Field of view (at 45 cm)	4.0x: 6.5 cm   5.0x: 6 cm   6.0x: 5.5 cm   7.0x: 5 cm
Depth of field (up to)	4.0x: 12.5 cm   5.0x: 11.5 cm   6.0x: 8.5 cm   7.0x: 7 cm
Mount system	TTL (Through the Lens)
Weight of two oculars with correction lenses	49 g
Loupe colour	Anthracite grey/Metallic grey
Loupe material	A Carbon enhanced composite housing and a titanium ring for changing magnification
Frames available	All