



# Kepler Reflekt

## Technical Specifications

Variant	Kepler Reflekt
Lenses	Plano or Prescription
BlueLight Protection (BLP)	Yes (optional add-on)
Prisma Lens System	Yes (optional add-on)
Fine-focus available	Yes
Magnification	4.0x   6.0x
Working distance	4.0x: From 44 to 62 cm   6.0x: From 46 to 68 cm
Field of view (at 45 cm)	4.0x: 7.2 cm   6.0x: 5.6 cm
Depth of field (up to)	4.0x: 13 cm   6.0x: 8.0 cm
Mount System	TTL (Through the lens)
Weight of two oculars with correction lenses	4.0x: 37.2 g   6.0x: 40.7 g
Loupe colour	Black housing with titanium front
Loupe material	Magnesium and titanium
Frames available	All