



# Total Pure & PowerGo

## Technical Specifications

Product name:	Total Pure v2 & PowerGo v2
Material	Aluminium
Weight	Total Intense: 33 g (Incl. Framepart & Cable) PowerGo: 155 g
Brightness settings	Low   Medium   High
Operating time	16h   12.5h   9h
Charging time	Approx. 3.5-5 hours
Light intensity	At 25 cm: 43,000 lux   At 35 cm: 21,500 lux
Light spot at working distance 40 cm	63 mm
Light colour	5,000 K
Colour rendering index	>90CRI
Cable type	Lamp Cable: Custom USB-C, 1.3 m Power Cord: USB-C 5W 2A
Battery type	21,700 Lithium-ion rechargeable battery
Battery rating	In: 5 VDC 2.0-2.5 A   Out: 3 VDC 0.5-2.0 A
Warranty	Total Pure: 3 years   PowerGo: 1 year
Product information	Flicker Free Light Rotatable belt clip Easy-to-operate Power button Over-current protection Over-voltage protection